


I

SYSTEM

PRODUCT CATALOG


IZEN
I M P L A N T



FIXTURE LINE UP


X-Narrow

Ø3.0
(mm)




R

Ø3.5
(mm)




R

Ø4.0
(mm)




R

Ø4.5
(mm)




W

Ø5.0
(mm)



W

Ø5.5
(mm)



• Special Options

INNER CONNECTION

Platform Switching

Prevents crestal bone loss after implant placement and provides esthetic result.

Two Conical contact points

Strong tightening force of dual contact between fixture and abutment results in less sink-down which leads to preventing screw loosening.



ZENEX PLUS I-System Fixture



Code	Article	Max. Diameter	Length
ZPSX3010	X-Narrow	Ø3.0	10.0mm
ZPSX3011	X-Narrow	Ø3.0	11.5mm
ZPSX3013	X-Narrow	Ø3.0	13.0mm
ZPSX3015	X-Narrow	Ø3.0	15.0mm



Code	Article	Max. Diameter	Length
ZPSM3508	Regular	Ø3.5	8.5mm
ZPSM3510	Regular	Ø3.5	10.0mm
ZPSM3511	Regular	Ø3.5	11.5mm
ZPSM3513	Regular	Ø3.5	13.0mm
ZPSM3515	Regular	Ø3.5	15.0mm



Code	Article	Max. Diameter	Length
ZPSR4007	Regular	Ø4.0	7.0mm
ZPSR4008	Regular	Ø4.0	8.5mm
ZPSR4010	Regular	Ø4.0	10.0mm
ZPSR4011	Regular	Ø4.0	11.5mm
ZPSR4013	Regular	Ø4.0	13.0mm
ZPSR4015	Regular	Ø4.0	15.0mm



Code	Article	Max. Diameter	Length
ZPSR4507	Regular	Ø4.5	7.0mm
ZPSR4508	Regular	Ø4.5	8.5mm
ZPSR4510	Regular	Ø4.5	10.0mm
ZPSR4511	Regular	Ø4.5	11.5mm
ZPSR4513	Regular	Ø4.5	13.0mm
ZPSR4515	Regular	Ø4.5	15.0mm



Code	Article	Max. Diameter	Length
ZPSW5007	Wide	Ø5.0	7.0mm
ZPSW5008	Wide	Ø5.0	8.5mm
ZPSW5010	Wide	Ø5.0	10.0mm
ZPSW5011	Wide	Ø5.0	11.5mm
ZPSW5013	Wide	Ø5.0	13.0mm
ZPSW5015	Wide	Ø5.0	15.0mm



Code	Article	Max. Diameter	Length
ZPSW5507	Wide	Ø5.5	7.0mm
ZPSW5508	Wide	Ø5.5	8.5mm
ZPSW5510	Wide	Ø5.5	10.0mm
ZPSW5511	Wide	Ø5.5	11.5mm
ZPSW5513	Wide	Ø5.5	13.0mm

1.2 Cover Screw



Code	Article	Height
CSSR40S	R/W	0.4mm
CSSR40M	R/W	1.4mm
CSSR40L	R/W	2.0mm

1.3 Healing Abutment



Code	Article	Diameter	Height
HASR402	M/R/W	Ø4.3	2.0mm
HASR403	M/R/W	Ø4.3	3.0mm
HASR404	M/R/W	Ø4.3	4.0mm
HASR405	M/R/W	Ø4.3	5.0mm
HASR452	M/R/W	Ø4.8	2.0mm
HASR453	M/R/W	Ø4.8	3.0mm
HASR454	M/R/W	Ø4.8	4.0mm
HASR455	M/R/W	Ø4.8	5.0mm
HASR522	M/R/W	Ø5.5	2.0mm
HASR523	M/R/W	Ø5.5	3.0mm
HASR524	M/R/W	Ø5.5	4.0mm
HASR572	M/R/W	Ø6.0	2.0mm
HASR573	M/R/W	Ø6.0	3.0mm
HASR574	M/R/W	Ø6.0	4.0mm

1.4 Temporary Abutment



Code	Article	Diameter	G/H	Connection
TASRH4510	M/R/W	Ø4.5	1.0mm	Hex
TASRN4510	M/R/W	Ø4.5	1.0mm	Non-Hex
TASRH4530	M/R/W	Ø4.5	3.0mm	Hex
TASRN4530	M/R/W	Ø4.5	3.0mm	Non-Hex



Code	Article	Diameter	G/H	P/H	Connection
CASRH45107	M/R/W	Ø4.5	1.0mm	7.0mm	Hex
CASRH45207	M/R/W	Ø4.5	2.0mm	7.0mm	Hex
CASRH45307	M/R/W	Ø4.5	3.0mm	7.0mm	Hex
CASRH45407	M/R/W	Ø4.5	4.0mm	7.0mm	Hex



Code	Article	Diameter	G/H	P/H	Connection
CASRN45107	M/R/W	Ø4.5	1.0mm	7.0mm	Non-Hex
CASRN45207	M/R/W	Ø4.5	2.0mm	7.0mm	Non-Hex
CASRN45307	M/R/W	Ø4.5	3.0mm	7.0mm	Non-Hex
CASRN45407	M/R/W	Ø4.5	4.0mm	7.0mm	Non-Hex



Code	Article	Diameter	G/H	P/H	Connection
CASRH52107	M/R/W	Ø5.2	1.0mm	7.0mm	Hex
CASRH52207	M/R/W	Ø5.2	2.0mm	7.0mm	Hex
CASRH52307	M/R/W	Ø5.2	3.0mm	7.0mm	Hex
CASRH52407	M/R/W	Ø5.2	4.0mm	7.0mm	Hex



Code	Article	Diameter	G/H	P/H	Connection
CASRN52107	M/R/W	Ø5.2	1.0mm	7.0mm	Non-Hex
CASRN52207	M/R/W	Ø5.2	2.0mm	7.0mm	Non-Hex
CASRN52307	M/R/W	Ø5.2	3.0mm	7.0mm	Non-Hex
CASRN52407	M/R/W	Ø5.2	4.0mm	7.0mm	Non-Hex



Code	Article	Diameter	G/H	P/H	Connection
CASRH57107	M/R/W	Ø5.7	1.0mm	7.0mm	Hex
CASRH57207	M/R/W	Ø5.7	2.0mm	7.0mm	Hex
CASRH57307	M/R/W	Ø5.7	3.0mm	7.0mm	Hex
CASRH57407	M/R/W	Ø5.7	4.0mm	7.0mm	Hex



Code	Article	Diameter	G/H	P/H	Connection
CASRN57107	M/R/W	Ø5.7	1.0mm	7.0mm	Non-Hex
CASRN57207	M/R/W	Ø5.7	2.0mm	7.0mm	Non-Hex
CASRN57307	M/R/W	Ø5.7	3.0mm	7.0mm	Non-Hex
CASRN57407	M/R/W	Ø5.7	4.0mm	7.0mm	Non-Hex

1.5 Angled Abutment



Code	Article	Diameter	G/H	Angle	Connection
AASRA45152	M/R/W	Ø4.5	2.0mm	15°	A-Type
AASRA45153	M/R/W	Ø4.5	3.0mm	15°	A-Type
AASRA45154	M/R/W	Ø4.5	4.0mm	15°	A-Type
AASRB45152	M/R/W	Ø4.5	2.0mm	15°	B-Type
AASRB45153	M/R/W	Ø4.5	3.0mm	15°	B-Type
AASRB45154	M/R/W	Ø4.5	4.0mm	15°	B-Type
AASRN45152	M/R/W	Ø4.5	2.0mm	15°	Non-Hex
AASRN45153	M/R/W	Ø4.5	3.0mm	15°	Non-Hex
AASRN45154	M/R/W	Ø4.5	4.0mm	15°	Non-Hex

1.6 FreeMilling Abutment



Code	Article	Diameter	G/H	Height	Connection
FMASRH45209	M/R/W	Ø4.5	2.0mm	9.0mm	Hex
FMASRN45209	M/R/W	Ø4.5	2.0mm	9.0mm	Non-Hex
FMASRH45409	M/R/W	Ø4.5	4.0mm	9.0mm	Hex
FMASRN45409	M/R/W	Ø4.5	4.0mm	9.0mm	Non-Hex



Code	Article	Diameter	G/H	Height	Connection
FMASRH57208	M/R/W	Ø5.7	2.0mm	8.0mm	Hex
FMASRN57208	M/R/W	Ø5.7	2.0mm	8.0mm	Non-Hex
FMASRH57408	M/R/W	Ø5.7	4.0mm	8.0mm	Hex
FMASRN57408	M/R/W	Ø5.7	4.0mm	8.0mm	Non-Hex



Code	Article	Diameter	G/H	Height	Connection
FMASRH52208	M/R/W	Ø5.2	2.0mm	8.0mm	Hex
FMASRN52208	M/R/W	Ø5.2	2.0mm	8.0mm	Non-Hex
FMASRH52408	M/R/W	Ø5.2	4.0mm	8.0mm	Hex
FMASRN52408	M/R/W	Ø5.2	4.0mm	8.0mm	Non-Hex

1.7 CCM Cast Abutment



Code	Article	Diameter	G/H	Connection
CCASRH4510	M/R/W	Ø4.5	1.0mm	Hex
CCASRN4510	M/R/W	Ø4.5	1.0mm	Non-Hex
CCASRH4530	M/R/W	Ø4.5	3.0mm	Hex
CCASRN4530	M/R/W	Ø4.5	3.0mm	Non-Hex


1.8 Components of Impression

Pick-Up Impression Coping

Listing Number D432520



Code	Article	Diameter	Length	Connection
PISR45S	M/R/W	Ø4.8	11.0mm	Hex
PISRN45S	M/R/W	Ø4.8	11.0mm	Non-Hex
PISR45L	M/R/W	Ø4.8	15.0mm	Hex
PISRN45L	M/R/W	Ø4.8	15.0mm	Non-Hex



Code	Article	Length
PIGPSRS	S Guide Pin (M/R/W)	11.0mm
PIGPSRL	L Guide Pin (M/R/W)	15.0mm



Code	Article	Diameter	Length	Connection
PISR40S	M/R/W	Ø4.0	11.0mm	Hex
PISRN40S	M/R/W	Ø4.0	11.0mm	Non-Hex
PISR40L	M/R/W	Ø4.0	15.0mm	Hex
PISRN40L	M/R/W	Ø4.0	15.0mm	Non-Hex

Transfer Impression Coping

Listing Number D432520



Code	Article	Diameter	Length	Connection
TISR45S	M/R/W	Ø4.8	9.0mm	Hex
TISRN45S	M/R/W	Ø4.8	9.0mm	Non-Hex
TISR45L	M/R/W	Ø4.8	13.0mm	Hex
TISRN45L	M/R/W	Ø4.8	13.0mm	Non-Hex



Code	Article	Length
TIGPSRS	S Guide Pin (M/R/W)	9.0mm
TIGPSRL	L Guide Pin (M/R/W)	13.0mm



Code	Article	Diameter	Length	Connection
TISR40S	M/R/W	Ø4.0	9.0mm	Hex
TISRN40S	M/R/W	Ø4.0	9.0mm	Non-Hex
TISR40L	M/R/W	Ø4.0	13.0mm	Hex
TISRN40L	M/R/W	Ø4.0	13.0mm	Non-Hex

Fixture Lab Analog



Code	Article	Diameter
FLASR400	R/W	Ø4.25

1.9

Multi Straight Abutment

Move into internal oral part by using exclusive Abutment Carrier (Code: MSACR48) [Listing Number D432520](#)

Code	Article	Diameter	G/H
MSASR4815	M/R/W	Ø4.8	1.5mm
MSASR4825	M/R/W	Ø4.8	2.5mm
MSASR4835	M/R/W	Ø4.8	3.5mm
MSASR4845	M/R/W	Ø4.8	4.5mm

1.10

Multi Angled Abutment

Connect by using exclusive Abutment Carrier (Code: MAACRMC) [Listing Number D432520](#)

Code	Article	Diameter	G/H	Angle	Connection
MAASRH481725	M/R/W	Ø4.8	2.5mm	17°	Hex
MAASRN481725	M/R/W	Ø4.8	2.5mm	17°	Non-Hex
MAASRH481735	M/R/W	Ø4.8	3.5mm	17°	Hex
MAASRN481735	M/R/W	Ø4.8	3.5mm	17°	Non-Hex
MAASRH481745	M/R/W	Ø4.8	4.5mm	17°	Hex
MAASRN481745	M/R/W	Ø4.8	4.5mm	17°	Non-Hex



Code	Article	Diameter	G/H	Angle	Connection
MAASRH483035	M/R/W	Ø4.8	3.5mm	30°	Hex
MAASRN483035	M/R/W	Ø4.8	3.5mm	30°	Non-Hex
MAASRH483045	M/R/W	Ø4.8	4.5mm	30°	Hex
MAASRN483045	M/R/W	Ø4.8	4.5mm	30°	Non-Hex

1.11

Components of Multi Abutment



Code	Product	Diameter
MTCSR48	Multi Temporary Cylinder	Ø4.8



Code	Product	Diameter
MPCCSR48	Multi Plastic Cast Cylinder	Ø4.8

[Listing Number D432520](#)

Code	Product	Diameter
MHCSR48	Multi Healing Cap	Ø4.8

[Listing Number D432520](#)

Code	Product	Diameter
MTICSR48	Multi Transfer Impression Coping	Ø4.8

[Listing Number D432520](#)

Code	Product
MSBSR4809	Multi ScanBody

[Listing Number D432520](#)

Code	Product	Diameter
MCCCSR48	Multi CCM Cast Cylinder	Ø4.8



Code	Product	Diameter
MTLCSR48	Multi Ti Link Cylinder	Ø4.8



Code	Product	Diameter
MPICSR48	Multi Pick-Up Impression Coping	Ø4.8

[Listing Number D432520](#)

Code	Product
MLASR48	Multi Lab Analog

[Listing Number D432520](#)

Code	Product
MSADSR20	Multi Straight Abutment Driver

[Listing Number D432520](#)

1.12

Ball Abutment

Ball Abutment



Code	Article	Diameter	G/H
BASR3510	M/R/W	Ø3.5	1.0mm
BASR3520	M/R/W	Ø3.5	2.0mm
BASR3530	M/R/W	Ø3.5	3.0mm
BASR3540	M/R/W	Ø3.5	4.0mm

1.13

ScanBody

Listing Number D432520



Code	Article	Diameter	Length
ISBR4310	R/W	Ø4.3	10.0mm

1.14

Ti Link Abutment



Code	Article	Diameter	G/H	P/H	Connection
TLASRH4510	M/R/W	Ø4.5	1.0mm	3.8mm	Hex
TLASRN4510	M/R/W	Ø4.5	1.0mm	3.8mm	Non-Hex
TLASRH4520	M/R/W	Ø4.5	2.0mm	3.8mm	Hex
TLASRN4520	M/R/W	Ø4.5	2.0mm	3.8mm	Non-Hex

Components of Ball Abutment



Code	Product
BASRCS	Ball Abutment Cap Set



Code	Product
BASRRS	Ball Abutment Retainer Set



Code	Product
BASROS	Ball Abutment O-ring Set

Listing Number D432520



Code	Product
BAALA	Ball Abutment Lab Analog

Listing Number D432520



Code	Product
BAD24	Ball Abutment Driver

Listing Number D432520



HQ & Factory
1, 2Dong, 26-32, Suworam 4-gil, Seotan-myeon,
Pyeongtaek-si, Gyeonggi-do, 17703
Republic of Korea

Overseas Business Department
10F 1dong #1011, #1012, #1013, #1014 Central-Tower II,
158, Magokseo-ro, Gangseo-gu, Seoul, 07788
Republic of Korea

Tel. +82-31-662-0657 (HQ)
+82-2-3664-3666 (Sales)
Fax +82-31-662-0658

izenimplant.com

Nouveta Healthcare Inc
R-45, 1st Floor Faraaz Complex
Ramesh Park, Laxmi Nagar
Delhi – 110092
GSTIN: 07AAC PQ1711A1ZF